

ZATA Standard

Name:

DoB:

Date:

Right Rx:.....

Pupil Diameter:OD..... OS.....

Record number

Age:

Time: 10:57:30 AM

Left Rx:.....

VA: OD OS

Threshold

30/24-2


Goldmann III white stimulus

Background 10 cd/m2



Blindspot Fixation Monitor

Central Fixation


RIGHT EYE



Threshold dB (0=3183cd/m2)



Total Defect dB



Pattern Defect dB

Presentations 54
Fovea Off



Fixation Losses 0/3

False Positives 0/5

False Negatives 0/3

Mean Defect 0.19 dB
Within normal limits

Std Dev Defect 0.59 dB
Within normal limits

Hemifield Test
Within normal limits

<5%
<2%
<1%
<0.5%

Henson 9000

Normative database Version 2.0

Software Version 4.0.114